

INFORMATION BROCHURE
FOR USERS AND RELATIVES

What is an electronic medicine dispenser?



Medicine dispensers

Do you find it difficult to remember to take your medications?

Then an electronic medicine dispenser can help you. It alerts you and doses out medicines at specific times.

There are different types of medicine dispensers. The one that is right for you will depend entirely on whether you need help taking your medications.

By using a medicine dispenser, you can take your medication at the right time without worrying about forgetting to take it, or taking the wrong ones at the wrong times. You also do not have to wait for a visit from the home care services and thus have greater freedom and flexibility in everyday life.

If you do not take your medication at the specified time, the Response Center will receive an alert. Then they will contact you to remind you to take your medicine. The Response Center can also contact the municipal health service if needed.



Automatic dosette

An automatic dosette ensures that you get the right medicine at the right time in an easy and safe way. The dosette alerts you when a dose of medicine should be taken.

An automatic dosette has 28 dosing compartments that alerts you by blinking and/or making a sound when the medicine is to be taken.

The dosette is durable and easy to bring with you wherever you want to go.

The dosette can be set to administer medicine one or more times a day.



Electronic medicine dispenser

Electronic medicine dispensers dispense medicines that are pre-dosed by the pharmacy. This ensures that you get the right medicine at the right time and notifies you when a dose of medicine should be taken. It alerts by a flashing light and making a sound.



To take the medicine, press the OK button and the medicine bag is fed out of the dispenser. The dispenser can be set to administer medicine once or several times a day.

The medicine dispenser shall be placed in a suitable place, with an outlet available. The medicine dispenser should not be placed on a sink or near other water sources, on an unstable surface (such as an unsteady table or chair), in direct sunlight or in the immediate vicinity of a heat source, such as a stove or fireplace.

If you are not at home when you need medicine, or if you are going to travel, you can pre-dispense the sachets from the dispenser. Travel Mode can usually only be activated by a caregiver/nursing staff, but you may be granted permission to use this.



Photo: Dignio

For more information, you can contact the health and care services in your municipality.



VIDEO – ELECTRONIC MEDICINE SUPPORT

You can scan this QR code with the camera on your smartphone or tablet to watch the film clip, “Electronic medicine support”



**RESPONSSENTERET (the Response Center)
– YOUR ALARM RECEIVING CENTER**

You can scan this QR code to see a presentation of Responscenteret (the Response Center)